



USERS GUIDE

Radiant QuikRAM

HDR and radioactive-materials management with built-in TG-43.

QuikRAM
RADIANT MEDICAL PHYSICS CONSULTING

InventorySource ModelsWell ChambersHDR Source ExchangesTG-43 CalcDecay TableReportsLOCAL-FIRST · SQLITE

Activate Radiant QuikRAM

Radiant QuikRAM requires a license key to run. Paste your **RADKEY~** key below to activate this install.

LICENSE KEY — the full key, with or without dashes

RADKEY-XXXXXX-XXXXXX-~

ACTIVATE QUIKRAM

Radiant QuikRAM verifies license keys entirely offline. No information leaves this server. Need a key? Contact Radiant Medical Physics Consulting.

Radiant QuikRAM — a tool of Radiant Medical Physics Consulting LLC. Runs 100% locally; no data leaves this machine. Source records, chamber-calibration values, and TG-43 calculations must be independently verified by a Qualified Medical Physicist before clinical use.

OVERVIEW

What QuikRAM is.

Radiant QuikRAM is HDR and radioactive-materials management for radiation oncology — source inventory with automatic isotope decay correction, the full HDR Source Exchange QA worksheet, well-chamber inventory with ADCL calibration values, a TG-43 line-source dose calculator, and printable date / AKS / activity decay tables — all in one local web app.

Who it's for

- HDR brachytherapy programs.
- Radiation safety officers managing sealed sources.
- Clinical physicists doing TG-43 reference checks.

QuikRAM
RADIANT MEDICAL PHYSICS CONSULTING

InventorySource ModelsWell ChambersHDR Source ExchangesTG-43 CalcDecay TableReportsLOCAL-FIRST · SQLITE

Activate Radiant QuikRAM

Radiant QuikRAM requires a license key to run. Paste your `RADKEY-...` key below to activate this install.

LICENSE KEY — the full key, with or without dashes

`RADKEY-XXXXXX-XXXXXX-...`

ACTIVATE QUIKRAM

Radiant QuikRAM verifies license keys entirely offline. No information leaves this server. Need a key? Contact Radiant Medical Physics Consulting.

Radiant QuikRAM — a tool of Radiant Medical Physics Consulting LLC. Runs 100% locally; no data leaves this machine.
Source records, chamber-calibration values, and TG-43 calculations must be independently verified by a Qualified Medical Physicist before clinical use.

SOURCE INVENTORY — DECAY-CORRECTED AKS TODAY.

Access and your source inventory.

Access

QuikRAM is deployed at quikram.radiant-mpc.com. Open the URL in any modern browser.

Add a source

1. From the dashboard, click *+ New source*.
2. Pick the isotope (Ir-192, Co-60, Cs-137, I-125, Pd-103, Y-90, etc.). QuikRAM uses the published half-life for decay corrections.
3. Enter the calibration date, the calibration air-kerma strength (or activity), and the certificate number.
4. Save. The source dashboard shows the AKS / activity decay-corrected to today.

QuikRAM
RADIANT MEDICAL PHYSICS CONSULTING

InventorySource ModelsWell ChambersHDR Source ExchangesTG-43 CalcDecay TableReportsLOCAL-FIRST · SQLITE

Activate Radiant QuikRAM

Radiant QuikRAM requires a license key to run. Paste your `RADKEY-...` key below to activate this install.

LICENSE KEY — the full key, with or without dashes

RADKEY-XXXXXX-XXXXXX-...

ACTIVATE QUIKRAM

Radiant QuikRAM verifies license keys entirely offline. No information leaves this server. Need a key? Contact Radiant Medical Physics Consulting.

Radiant QuikRAM — a tool of Radiant Medical Physics Consulting LLC. Runs 100% locally; no data leaves this machine.
Source records, chamber-calibration values, and TG-43 calculations must be independently verified by a Qualified Medical Physicist before clinical use.

SOURCE DETAIL — DECAY CORRECTION TO TODAY.

HDR Source Exchange QA.

The Source Exchange workflow walks you through every section of the quarterly QA in order, with auto-computed metrics where the math is mechanical.

1. **ODH acceptance checks** on the new source.
2. **Chamber constancy** on the prior source (before removal).
3. **New source calibration** — well-chamber reading, correction factors, calculated AKS, comparison to vendor cert.
4. **Timer accuracy, reproducibility, and linearity** with an auto-computed R-squared plot.
5. **Housing survey** and **source-positioning film** upload.

At the end, click *Generate report* for a complete QA PDF ready for the source binder.

QuikRAM
RADIANT MEDICAL PHYSICS CONSULTING

InventorySource ModelsWell ChambersHDR Source ExchangesTG-43 CalcDecay TableReportsLOCAL-FIRST · SQLITE

Activate Radiant QuikRAM

Radiant QuikRAM requires a license key to run. Paste your `RADKEY-...` key below to activate this install.

LICENSE KEY — the full key, with or without dashes

RADKEY-XXXXXX-XXXXXX-...

ACTIVATE QUIKRAM

Radiant QuikRAM verifies license keys entirely offline. No information leaves this server. Need a key? Contact Radiant Medical Physics Consulting.

Radiant QuikRAM — a tool of Radiant Medical Physics Consulting LLC. Runs 100% locally; no data leaves this machine.
Source records, chamber-calibration values, and TG-43 calculations must be independently verified by a Qualified Medical Physicist before clinical use.

HDR SOURCE EXCHANGE QA — SECTIONED WORKSHEET.

Well-chamber inventory + TG-43 calculator.

Well-chamber inventory

Each well chamber is a tracked asset with manufacturer, serial number, ADCL calibration coefficient, calibration date, and the cross-calibration history. The dashboard shows when a chamber's ADCL is due.

TG-43 line-source calculator

An independent dose-rate reference check. Pick the source, the distance, and the off-axis angle; QuikRAM applies the published $g(r)$, $F(r, \theta)$, and anisotropy data and returns the dose-rate at point. Useful as a second check on planning-system output, especially during commissioning.

QuikRAM
RADIANT MEDICAL PHYSICS CONSULTING

InventorySource ModelsWell ChambersHDR Source ExchangesTG-43 CalcDecay TableReportsLOCAL-FIRST · SQLITE

Activate Radiant QuikRAM

Radiant QuikRAM requires a license key to run. Paste your `RADKEY-...` key below to activate this install.

LICENSE KEY — the full key, with or without dashes

RADKEY-XXXXXX-XXXXXX-...

ACTIVATE QUIKRAM

Radiant QuikRAM verifies license keys entirely offline. No information leaves this server. Need a key? Contact Radiant Medical Physics Consulting.

Radiant QuikRAM — a tool of Radiant Medical Physics Consulting LLC. Runs 100% locally; no data leaves this machine.
Source records, chamber-calibration values, and TG-43 calculations must be independently verified by a Qualified Medical Physicist before clinical use.

WELL-CHAMBER DETAIL — ADCL COEFFICIENT AND HISTORY.

Decay tables and binder reports.

Date / AKS / activity decay tables

Pick a source and a horizon (30, 60, 90, 180 days) and QuikRAM prints a day-by-day decay table — the table the source binder needs and the table that backs the day's planning calculations.

Source-exchange report

The exchange workflow rolls every section into a Radiant-branded PDF: acceptance checks, constancy, calibration, timer / reproducibility / linearity, housing survey, and film image, with the reviewer signature block on the last page.

QuikRAM
RADIANT MEDICAL PHYSICS CONSULTING

InventorySource ModelsWell ChambersHDR Source ExchangesTG-43 CalcDecay TableReportsLOCAL-FIRST · SQLITE

Activate Radiant QuikRAM

Radiant QuikRAM requires a license key to run. Paste your `RADKEY-...` key below to activate this install.

LICENSE KEY — the full key, with or without dashes

RADKEY-XXXXXX-XXXXXX-...

ACTIVATE QUIKRAM

Radiant QuikRAM verifies license keys entirely offline. No information leaves this server. Need a key? Contact Radiant Medical Physics Consulting.

Radiant QuikRAM — a tool of Radiant Medical Physics Consulting LLC. Runs 100% locally; no data leaves this machine. Source records, chamber-calibration values, and TG-43 calculations must be independently verified by a Qualified Medical Physicist before clinical use.

DECAY TABLE — DATE, AKS, ACTIVITY COLUMNS.

License, safety, and how to reach us.

License terms

QuikRAM is licensed per HDR program. The deployed instance at quikram.radiant-mpc.com is provided as a hosted evaluation; a self-hosted on-prem build is available on request.

Evaluation status

QuikRAM is an unvalidated analysis aid. Every TG-43 calculation and every source-exchange QA conclusion must be independently reviewed by a Qualified Medical Physicist before clinical use.

Contact

Live URL: <https://quikram.radiant-mpc.com/>

Email: info@Radiant-MPC.com

Website: radiant-mpc.com

Evaluation / decision-support only. Every Radiant tool is an aid to a Qualified Medical Physicist. All outputs must be independently reviewed before any clinical use. Do not transmit PHI through evaluation builds.