



Radiant QuikShare

HIPAA-conscious, AES-256 DICOM transfer with audit trail.

IN DEVELOPMENT

QuikShare is on the build queue. The DICOM scan + multi-patient guard, .7z header-encrypted archive layer, OneDrive delivery, and HIPAA accounting-of-disclosures log are scoped; deployment targets a workstation-local app with OneDrive for Business as the secured transfer surface.

WHAT IT DOES

Package a complete planning dataset (CT, RT Plan, RT Dose, RT Structure Set), encrypt it with AES-256 (header encryption, so even the filenames inside are hidden), generate a 20-character random passphrase, upload the archive to OneDrive for Business, create an expiring sharing link, and log the disclosure for HIPAA accounting. The multi-patient guard refuses to proceed if files from multiple patients are mixed in.

KEY CAPABILITIES

Patient-Safe Scanning

DICOM inventory; refuses on mixed-patient input.

AES-256 Encryption

.7z archive with header encryption (filenames hidden).

Random Passphrase

20-character cryptographic passphrase with verbal formatting.

OneDrive Delivery

Upload + expiring sharing link.

Integrity Verification

SHA-256 hash; drafts the recipient email.

HIPAA Audit Trail

Dual-location disclosure log; passphrases stored as salted hashes.

WHO IT'S FOR

- Clinical medical physicists
- Inter-center plan transfers

- Second-opinion and consult workflows
- HIPAA accounting-of-disclosures

HOW TO GET ACCESS

QuikShare is not yet deployed. To request a pilot or join the launch list, email info@Radiant-MPC.com.

SUPPORT

Questions, feature requests, or sales inquiries: info@Radiant-MPC.com

Web: radiant-mpc.com