



USERS GUIDE

Radiant QuikLog

Equipment master and calibration tracker for the QA program.

QuikLog
RADIANT MEDICAL PHYSICS CONSULTING

Inventory + Add Asset Inventory Report Due Soon LOCAL-FIRST - SQLite

Equipment Inventory

12 assets tracked — May 23, 2026. Items overdue or due within 60 days are flagged below.

6
OK

1
DUE SOON

4
OVERDUE

1
NO DATE

+ Add Asset

Printable Inventory

Due in Next 60 Days

Status	Type	Asset	S/N	Location	Last Cal	Cal Due	Days
Overdue	Source-Strength Meter	Standard Imaging SourceCheck	SC-1144	HDR Vault	2025-03-29	2026-03-29	-55
Overdue	Ion Chamber	PTW Farmer 30013	F-9981	Physics Office	2025-04-18	2026-04-18	-35
Overdue	Well Chamber	Standard Imaging HDR1000+	A200788	HDR Vault	2025-05-08	2026-05-08	-15
Overdue	CT Simulator	Siemens SOMATOM go.Sim	SS-2024-019	CT Sim	2025-06-07	2026-05-18	-5
Due Soon	Linear Accelerator	Elekta Versa HD	VHD-4112	Vault 2	2025-07-27	2026-06-07	15
OK	HDR Unit	Elekta Flexitron	FX-2210	HDR Vault	2025-11-24	2026-11-24	185
OK	Survey Meter / Geiger	Ludlum 9DP-1	LD-3382	HDR Vault	2025-11-19	2026-11-24	185
OK	Linear Accelerator	Varian TrueBeam	TB-7501	Vault 1	2026-05-03	2027-05-03	345
OK	Electrometer	Standard Imaging SUPERMAX	X212455	Physics Office	2025-10-25	2027-10-25	520
OK	Ion Chamber	Exradin A12	A12-2207	Physics Office	2026-01-31	2028-01-31	610

OVERVIEW

What QuikLog is.

Radiant QuikLog is an equipment master and calibration tracker for radiation oncology — linacs, CT simulators, HDR units, ion chambers, electrometers, well chambers, and survey meters, in one local inventory sorted by what's due first.

The single source of truth

Every piece of QA-relevant equipment lives in one searchable list, with serial numbers, manufacturers, install dates, and the full calibration history per asset. Color-coded due-soon and overdue badges put attention exactly where it's needed.

Who it's for

- Clinical physicists and QA programs.
- Equipment masters for state inspections and J-CAHO visits.
- ADCL coordination and chamber-bank management.

QuikLog

RADIANT MEDICAL PHYSICS CONSULTING

Inventory+ Add AssetInventory ReportDue SoonLOCAL-FIRST · SQLITE

Equipment Inventory

12 assets tracked — May 23, 2026. Items overdue or due within 60 days are flagged below.

6OK

1DUE SOON

4OVERDUE

1NO DATE

+ Add Asset

Printable Inventory

Due in Next 60 Days

Status	Type	Asset	S/N	Location	Last Cal	Cal Due	Days
Overdue	Source-Strength Meter	Standard Imaging SourceCheck	SC-1144	HDR Vault	2025-03-29	2026-03-29	-55
Overdue	Ion Chamber	PTW Farmer 30013	F-9981	Physics Office	2025-04-18	2026-04-18	-35
Overdue	Well Chamber	Standard Imaging HDR1000+	A200788	HDR Vault	2025-05-08	2026-05-08	-15
Overdue	CT Simulator	Siemens SOMATOM go.Sim	SS-2024-019	CT Sim	2025-06-07	2026-05-18	-5
Due Soon	Linear Accelerator	Elekta Versa HD	VHD-4112	Vault 2	2025-07-27	2026-06-07	15
OK	HDR Unit	Elekta Flexitron	FX-2210	HDR Vault	2025-11-24	2026-11-24	185
OK	Survey Meter / Geiger	Ludlum 9DP-1	LD-3382	HDR Vault	2025-11-19	2026-11-24	185
OK	Linear Accelerator	Varian TrueBeam	TB-7501	Vault 1	2026-05-03	2027-05-03	345
OK	Electrometer	Standard Imaging SUPERMAX	X212455	Physics Office	2025-10-25	2027-10-25	520
OK	Ion Chamber	Exradin A12	A12-2207	Physics Office	2026-01-23	2028-01-23	610

ASSET LIST — DUE-SOON AND OVERDUE FRONT-AND-CENTER.

Open QuikLog, add your first asset.

Access

QuikLog is deployed at quiklog.radiant-mpc.com. Open the URL in any modern browser. There is nothing to install on the workstation.

Sign-in

QuikLog uses a per-instance access password supplied with your agreement. Enter the password once and the session persists in the browser.

Add an asset

1. From the dashboard, click *+ New Asset*.
2. Set the asset type — Linac, CT Sim, HDR, Ion Chamber, Electrometer, Well Chamber, Survey Meter, or Other.
3. Fill in manufacturer, model, serial number, location, install date, and the recurring calibration interval.
4. Save. The asset appears on the dashboard with its first next-due date.

QuikLog
RADIANT MEDICAL PHYSICS CONSULTING

Inventory + Add Asset Inventory Report Due Soon LOCAL-FIRST · SQLITE

Inventory / Varian TrueBeam

Varian TrueBeam
Linear Accelerator · S/N TB-7501 · Vault 1

OK · 345 d

Details

Type	Linear Accelerator
Manufacturer / Model	Varian / TrueBeam
Serial number	TB-7501
Location	Vault 1
Status	In Service
Date of first use	2019-06-01
Last cal date	2026-05-03
Cal due date	2027-05-03
Cal interval (mo)	12
Notes	Primary photon/electron linac, IMRT/VMAT/SRS.

Edit asset+ Record calibrationDelete

Calibration history

ASSET DETAIL — CALIBRATION HISTORY AND NEXT-DUE.

Log a calibration event.

The calibration event is the unit of work in QuikLog.

1. Open the asset, click *+ Log calibration*.
2. Enter the calibration date, the next-due date (defaults from interval), the performer, and the ADCL or vendor certificate number.
3. Attach the certificate PDF if you have one — it travels with the event in subsequent reports.
4. Save. The event joins the asset's calibration ledger, and the asset's next-due date updates automatically.

Audit trail. Calibration events are append-only — corrections are made by logging a corrective event, never by editing history.

QuikLog
RADIANT MEDICAL PHYSICS CONSULTING

Inventory + Add Asset Inventory Report Due Soon LOCAL-FIRST - SQLite

Inventory / Varian TrueBeam / Record calibration

Record calibration

Linear Accelerator - Varian TrueBeam

Cal date

Next due date auto-filled from interval; override as needed

Performed by

Cert / ADCL #

Notes (optional)

Radiant QuikLog — a tool of Radiant Medical Physics Consulting LLC. Runs 100% locally; no data leaves this machine.
Calibration warning thresholds must be reviewed and approved by a Qualified Medical Physicist or Radiation Safety Officer.

NEW CALIBRATION EVENT — DATE, PERFORMER, ADCL NUMBER, CERTIFICATE UPLOAD.

The due-soon dashboard.

The due-soon report is the daily-glance view: every asset whose next-due date falls in the next N days, sorted soonest-first. The default window is 60 days; adjust with the *days* filter.

- Red rows are **overdue**.
- Amber rows are **due within 14 days**.
- Yellow rows are **due within 30 days**.

The dashboard is the right view for the start of a QA week — everything that needs attention is on one screen, with a one-click hop into the asset to log the cal once it's done.

QuikLog

RADIANT MEDICAL PHYSICS CONSULTING

Inventory+ Add AssetInventory ReportDue SoonLOCAL-FIRST · SQLITE

Inventory / Report

Calibrations Due in Next 60 Days

Print / Save as PDF

Radiant QuikLog · generated May 23, 2026 · 5 entries

Status	Type	Manufacturer / Model	S/N	Location	Last Cal	Cal Due	Days
Overdue	Source-Strength Meter	Standard Imaging SourceCheck	SC-1144	HDR Vault	2025-03-29	2026-03-29	-55
Overdue	Ion Chamber	PTW Farmer 30013	F-9981	Physics Office	2025-04-18	2026-04-18	-35
Overdue	Well Chamber	Standard Imaging HDR1000+	A200788	HDR Vault	2025-05-08	2026-05-08	-15
Overdue	CT Simulator	Siemens SOMATOM go.Sim	SS-2024-019	CT Sim	2025-06-07	2026-05-18	-5
Due Soon	Linear Accelerator	Elekta Versa HD	VHD-4112	Vault 2	2025-07-27	2026-06-07	15

Print this page or use your browser's Save as PDF to capture the report. Calibration warning thresholds must be reviewed and approved by a Qualified Medical Physicist or Radiation Safety Officer.

Radiant QuikLog — a tool of Radiant Medical Physics Consulting LLC. Runs 100% locally; no data leaves this machine. Calibration warning thresholds must be reviewed and approved by a Qualified Medical Physicist or Radiation Safety Officer.

DUE-SOON REPORT — COLOR-CODED BY URGENCY WINDOW.

Printable inventory + due-soon PDFs.

Two reports cover the binder side of the workflow.

Full inventory PDF

The complete asset master with serial numbers, location, calibration status, and next-due — formatted for the QA binder and for handing to state inspectors.

Due-soon PDF

The filtered list for the next N days, ready for the cal-week huddle.

Both reports inherit your practice name and any logo you uploaded in *Settings*.

QuikLog

RADIANT MEDICAL PHYSICS CONSULTING

Inventory + Add Asset Inventory Report Due Soon LOCAL-FIRST - SQLITE

Inventory / Report

Equipment Inventory

Radiant QuikLog · generated May 23, 2026 · 12 entries

Print / Save as PDF

Status	Type	Manufacturer / Model	S/N	Location	Last Cal	Cal Due	Days
Overdue	Source-Strength Meter	Standard Imaging SourceCheck	SC-1144	HDR Vault	2025-03-29	2026-03-29	-55
Overdue	Ion Chamber	PTW Farmer 30013	F-9981	Physics Office	2025-04-18	2026-04-18	-35
Overdue	Well Chamber	Standard Imaging HDR1000+	A200788	HDR Vault	2025-05-08	2026-05-08	-15
Overdue	CT Simulator	Siemens SOMATOM go.Sim	SS-2024-019	CT Sim	2025-06-07	2026-05-18	-5
Due Soon	Linear Accelerator	Elekta Versa HD	VHD-4112	Vault 2	2025-07-27	2026-06-07	15
OK	HDR Unit	Elekta Flexitron	FX-2210	HDR Vault	2025-11-24	2026-11-24	185
OK	Survey Meter / Geiger	Ludlum 9DP-1	LD-3382	HDR Vault	2025-11-19	2026-11-24	185
OK	Linear Accelerator	Varian TrueBeam	TB-7501	Vault 1	2026-05-03	2027-05-03	345
OK	Electrometer	Standard Imaging SUPERMAX	X212455	Physics Office	2025-10-25	2027-10-25	520
OK	Ion Chamber	Exradin A12	A12-2207	Physics Office	2026-01-23	2028-01-23	610
OK	Treatment Planning Workstation	Varian Eclipse 16.1 workstation	EW-22	TPS room	2026-02-12	2028-02-12	630
No Date	Survey Meter / Geiger	Fluke 451P	451-008	Physics Office	—	—	—

Print this page or use your browser's Save as PDF to capture the report. Calibration warning thresholds must be reviewed and approved by a Qualified Medical Physicist or Radiation Safety Officer.

INVENTORY REPORT — PRINTABLE PDF.

License, safety, and how to reach us.

License terms

QuikLog is licensed per practice. The deployed instance at quiklog.radiant-mpc.com is provided as a hosted evaluation; a self-hosted on-prem build is available on request.

Privacy

QuikLog stores equipment master data only — no patient identifiers and no PHI. Calibration certificate attachments are stored in the application's data directory.

Contact

Live URL: <https://quiklog.radiant-mpc.com/>

Email: info@Radiant-MPC.com

Website: radiant-mpc.com

Evaluation / decision-support only. Every Radiant tool is an aid to a Qualified Medical Physicist. All outputs must be independently reviewed before any clinical use. Do not transmit PHI through evaluation builds.