

DESIGN SPEC • QUIKBOLUS PHASE 2.5 • DRAFT V0.1

QuikBolus Phase 2.5 — Mold Stretch + Profile Expansion + Engineering Followups

Spec version: 0.1 (draft) **Date:** 2026-05-24 **Author:** Radiant Medical Physics Consulting / Claude Opus 4.7
Status: Design draft — ready for scoping; not yet committed to a release window. **Parent spec:** [design-docs/quikbolus-phase-2-spec.md](#) (defines Phases 2.1–2.5; this doc supersedes that doc's Phase 2.5 section with what we actually learned shipping 2.3 + 2.4) **Scope:** Three feature tracks + four engineering followups surfaced by the 2.3/2.4 production deploy. Direct-print profile work originally planned for Phases 2.1/2.2 remains pending and is rolled into the 2.5 scope here.

1. Executive summary

Phases 2.3 and 2.4 shipped the mold-design workflow end-to-end: [mold_design.py](#) engine, Flask routes, Three.js UI, registration pegs, draft angle, identification embossing, undercut detection, validation panel, PDF cast-info sheet, 3MF export with a Bambu Studio profile. All live at <https://quikbolus.radiant-mpc.com/mold> (commit [16df2de](#)).

Phase 2.5 takes the remaining work in two shapes:

1. **Mold stretch features that 2.4 explicitly deferred:** stepped parting planes for complex undercut anatomy.
 2. **Profile system expansion:** multi-printer 3MF profiles (Orca, PrusaSlicer, Cura — beyond the Bambu profile shipped in 2.4) for BOTH the mold workflow and the direct-print workflow; per-clinic admin overrides for default profile choices; and the printer-picker UI that was originally Phase 2.1.
 3. **Engineering followups surfaced by 2.4 in production:** a real wall-thickness check (placeholder today), undercut marker clustering (per-face today; noisy on dense overhangs), memory profile + targeted optimization (we OOM'd at 2 GB on real-anatomy mold builds before bumping Render to Pro 4 GB).
-

2. Motivation — what 2.4 production use surfaced

ISSUE FROM 2.3/2.4 LIVE USE	WHAT 2.5 DOES ABOUT IT
Single flat parting plane can't capture jaw + shoulder undercuts on H&N anatomy	Stepped parting plane lets the split follow anatomy at multiple Z levels
Wall-thickness check shows "not checked" in the validation panel	Real ray-cast sampler reports min wall thickness; warns below 2 mm
Undercut overlay renders one red sphere per offending face — chest-wall meshes can show 100+ markers, drowning the operator	Cluster face centroids into regions; one labelled marker per cluster
Render Standard tier (2 GB) OOM'd on real chest-wall mold builds; required upgrade to Pro (4 GB) at higher cost	Memory profile the engine; replace per-pixel text-mesh union with single shapely-polygon extrusion; chunk the cavity boolean for huge meshes
Only Bambu Studio gets a baked 3MF profile; clinics on Prusa MK4, Raise3D, Ultimaker, or OrcaSlicer get a generic STL	Per-printer profile bake — Orca, PrusaSlicer, Cura — same family of presets
No way for a multi-site clinic to set their own default print profile for everyone	Admin overrides: clinic-level defaults editable through Client Setup
Direct-print workflow has no printer picker at all; operator must know the right settings	Direct-print printer picker (was Phase 2.1 in the original spec)

3. Scope

In scope

- **3.1** Stepped parting planes for the mold workflow
- **3.2** Multi-printer 3MF profile bake for mold AND direct-print outputs (Bambu / Orca / PrusaSlicer / Cura)
- **3.3** Per-clinic admin overrides for default print profiles
- **3.4** Direct-print printer picker UI (the work originally in Phase 2.1)
- **3.5** Real wall-thickness sampler in the validation panel
- **3.6** Undercut marker clustering
- **3.7** Memory optimization — text-mesh single-extrusion path + boolean chunking for large cavities
- **3.8** Engineering hygiene: header nav (`/mold` reachable from inside the design + direct-print flows); session-level diagnostics endpoint (peak RAM, boolean engine times) for support cases

Out of scope (Phase 3.0+)

- Direct LAN printing (QuikBolos → Bambu MQTT / Klipper / OctoPrint → printer) — separate spec; networking + credentials + status polling
- Auto-orient bolus to anatomy from DICOM landmarks — separate spec

- HDR cap bolus templates / parametric anatomy library — Phase 3 product
- Tissue compensator / wedge generation — distinct product
- Multi-material bolus (rigid backing + flex contact surface) — research

4. Stepped parting planes (3.1)

A single horizontal Z-plane can't reliably split anatomy with multiple undercut levels — jaw + shoulder + ear is a classic case. Stepped parting "plane" is actually a staircase of Z heights connected by short vertical risers, traced through the bolus so each step lies above one cluster of undercuts and below the next.

4.1 Variables

VARIABLE	CHOICES
Step mode	Single plane (default — same as 2.3/2.4) · Stepped
Step count	2 (default — front and back halves) · 3 · 4+
Step positions	Auto-suggested from undercut detection · Manual click-to-place per step
Step transition style	Vertical wall (default — short Z riser between steps) · Ramped (sloped between steps; harder to print but cleaner cast surface)
Step minimum width	Default 8 mm (steps narrower than this collapse to a continuous slope)

4.2 Algorithm sketch

1. Run undercut detection at the current single Z choice; cluster offending faces into regions in XY (k-means or DBSCAN with eps ≈ 20 mm)
2. For each cluster, compute the "preferred Z" that minimizes its undercuts in isolation
3. Connect those Z values into a step staircase ordered along the bolus's principal XY axis
4. Replace the single `parting_z` cutting plane with a swept surface — each region uses its preferred Z, with vertical risers between
5. Boolean-split the mold using this swept surface instead of `slice_mesh_plane`. trimesh doesn't ship this directly; needs a custom `slice_with_swept_plane()` built from manifold3d boolean ops

4.3 Validation additions

- Each step ≥ `step_minimum_width` (else collapse / warn)
- No vertical riser exceeds 20 mm (would create unprintable overhang on the mold itself)
- Overall step count doesn't exceed `step_count` (warn if undercut clustering wants more steps than the operator allowed)

4.4 UI

- Step Mode dropdown in the parameters panel

- When "Stepped" selected, a new sub-panel appears with step-count and per-step Z entries
- 3D preview overlays the swept parting surface in semi-transparent yellow so the operator can see the staircase before committing

5. Multi-printer 3MF profiles (3.2)

Phase 2.4 shipped a Bambu Studio profile JSON embedded in the 3MF. Phase 2.5 adds three more profile families so clinics on non-Bambu printers get the same one-click experience.

5.1 Printer families to support

FAMILY	PROFILE SCHEMA
Bambu Studio (already shipped 2.4)	<code>Metadata/project_settings.config</code> JSON
OrcaSlicer	<code>Metadata/project_settings.config</code> JSON (similar schema; auto-detect by <code>profile_id</code> field)
PrusaSlicer	<code>Metadata/Slic3r_PE.config</code> INI-style key=value
Cura	<code>Metadata/Cura.proj</code> JSON

5.2 Per-printer profile content

SETTING	MOLD (RIGID)	DIRECT-PRINT BOLUS (TPU 95A)
Material	PETG (default) or PLA	TPU 95A
Layer height	0.16 mm	0.20 mm
Wall loops	4	3
Top/bottom shells	4	4
Infill density	35% gyroid	100% rectilinear
Nozzle temp	230-240 °C	240 °C
Bed temp	60-70 °C	55 °C
Supports	Normal at sprue/vent overhangs	Tree-organic at rim
Print speed (outer wall)	35 mm/s	25 mm/s

These match what the QuikBolus Print Operator Guide (the Bambu P2S PDF in `app/users-guides/quikbolus-print-guide.pdf`) recommends.

5.3 UI

The existing export dropdown grows:

```
Export format:
Two STLs (generic)
3MF -- Bambu Studio (P1S / P1P / X1C / X1E / A1 / H2D)
3MF -- OrcaSlicer
3MF -- PrusaSlicer (MK4 / MK4S / MK3S+ / XL / Mini+)
3MF -- Cura (Ultimaker S3 / S5 / S7)
```

Picking a slicer also reveals a "Target printer" sub-dropdown for the specific model within that ecosystem (drives build-volume fit check).

5.4 Schema durability

Slicer config schemas evolve. Phase 2.5 ships:

- Tested-against versions in `app/printer_profiles/` : Bambu Studio 1.9, OrcaSlicer 2.1, PrusaSlicer 2.7, Cura 5.6
- A `last_tested` field embedded in each profile so the UI can show "Last tested: Bambu Studio 1.9" — operator knows when to expect drift
- Profiles structured as a thin JSON template (mold + direct-print variants per printer) that's straightforward to update independent of the engine code

6. Per-clinic admin overrides (3.3)

A multi-clinic operator (running QuikBolus for multiple sites) shouldn't have to re-pick the default printer + material + sprue size for every case. Phase 2.5 makes the printer profile + mold defaults editable at the clinic level.

6.1 Variables surfaced

SETTING	PER-SITE DEFAULT
Default casting material	Flexibolus / Ecoflex 30 / Dragon Skin 10 / Custom
Default printer family	Bambu / Orca / Prusa / Cura
Default printer model within family	P1S / X1C / MK4 / etc.
Default mold parameters	wall, floor, headroom, sprue count, vent count, peg config
Default direct-print parameters	infill, layer height, walls/shells (for non-clinical or quick proofs)
Default identification embossing	enabled? format string template?

6.2 Storage

Per-clinic JSON record in the existing Client Setup table (`clients.html` in the marketing-site/app admin). New `default_quikbolus_settings` blob column on the clients row — Worker passes it down to QuikBolus on authenticated sessions.

6.3 UI surface

- New tab in `app/admin/clients.html` (the existing Client Setup UI) per client: **Default Settings** → **QuikBolos** → **Edit**
- When a customer signs into QuikBolos, their session pulls these defaults; the `/mold` and `/` (direct-print) param forms initialize to the per-clinic values instead of the hard-coded ones in `MoldParams`
- Per-case override still works — the per-clinic defaults are just the starting point

6.4 Why this matters

The single-clinic case is fine with hard-coded defaults; once you have 2+ clinics with different printer fleets the friction is real (every case starts with "wait, why is it set to Bambu when we run Prusa here?"). This is the lowest-effort feature in 2.5 with the highest UX payoff for multi-site customers.

7. Direct-print printer picker (3.4)

Originally Phase 2.1 in the parent spec; never shipped. Now folds into 2.5 alongside the multi-printer profile work.

7.1 Scope

- Target Printer dropdown on the main `/` upload page (direct-print workflow), parallel to the mold workflow's picker
- Build-volume fit check on the bolus mesh against the chosen printer
- Recommended-settings panel showing the values from §5.2 (direct-print column) for the picked printer
- 3MF export option with the picked printer's profile baked in
- STL export still available for "generic / no-slicer-preference"

7.2 UI

Mirror the mold-design page's printer-pick UX so operators see the same shape regardless of which workflow they're in.

8. Real wall-thickness check (3.5)

2.4's validation panel shows "Wall thickness: not checked" because trimesh 4.12 doesn't ship a sampler. 2.5 implements one.

8.1 Approach

For each cavity-facing face on the mold's interior surface, cast a ray along the inward normal (into the mold material) and measure distance to the nearest exterior face. Report:

- Min wall thickness across all samples
- Median wall thickness

- Histogram bins (< 2 mm / 2-3 mm / 3-5 mm / > 5 mm)

8.2 Validation rules

MIN WALL OBSERVED	SEVERITY
≥ 3.0 mm	OK
2.0 - 3.0 mm	Warning ("thin region — verify pour pressure won't blow it out")
< 2.0 mm	Error ("wall too thin — increase <code>wall_thickness_mm</code> or simplify cavity geometry")

8.3 Performance

Ray cast for every interior face is expensive — chest-wall molds can have 100k+ faces. Subsample: pick a uniform grid of ~5000 points on the cavity surface; cast rays from those only. Statistically representative, runs in < 2 seconds.

9. Undercut clustering (3.6)

2.4 ships one red sphere per offending face. On real anatomy that's visually noisy. 2.5 clusters them.

9.1 Approach

1. Run undercut detection (returns per-face centroids + severity)
2. DBSCAN clustering on XY centroids (eps = 15 mm, min_samples = 3)
3. For each cluster, compute the bounding box and severity (max across member faces)
4. Render one marker per cluster: a small flag at the cluster centroid, colored by severity, with a count label "23 faces"

9.2 UI

- Replace per-face red spheres with one marker per cluster
- Click a cluster marker → opens a side panel listing the member faces with their individual severity
- "Show all faces" toggle to revert to per-face mode for the debug case

10. Memory optimization (3.7)

We OOM'd at 2 GB on real-anatomy mold builds before bumping Render to 4 GB Pro. Two targeted wins:

10.1 Text-mesh path

Today: PIL rasterizes glyphs → per-pixel box meshes → manifold union. Cost: 200-500 boxes per glyph × 5-10 glyphs = thousands of meshes unioned together. Memory peak: ~400 MB.

Fix: rasterize → trace pixel boundaries to a shapely polygon → single 3D extrusion via

`trimesh.creation.extrude_polygon`. One mesh instead of thousands. Memory peak: ~10 MB. Output is

identical.

10.2 Boolean chunking for huge cavities

Today: full anatomy mesh → single boolean subtract against the block. Cost: manifold3d allocates ~5× input mesh as scratch space; a 100 MB patient mesh peaks at 500 MB during the boolean.

Fix: split the bolus into octree chunks (typically 8 chunks for chest-wall), boolean-subtract each chunk in sequence, union the results. Peak memory $\approx$ chunk-size $\times$ 5. For chest-wall this brings peak from ~800 MB to ~200 MB.

10.3 Telemetry endpoint

New `GET /mold/api/diag` returns per-step peak RAM + boolean engine times for the most recent build. Helps diagnose "this case OOMed" tickets without instrumenting Render itself.

11. UX flow updates

11.1 Mold workflow

```
DICOM RT → bolus mesh (existing) → mold-design page
├─ Casting material
├─ Step mode (single | stepped) ← NEW
├─ Parting plane(s) – interactive ← stepped variant ← NEW
├─ Wall / floor / headroom (existing)
├─ Sprue / vent (existing)
├─ Registration pegs (existing 2.4)
├─ Draft / embossing (existing 2.4)
├─ Export format
│   ├─ Two STLs
│   ├─ 3MF – Bambu Studio (existing 2.4)
│   ├─ 3MF – Orca ← NEW
│   ├─ 3MF – PrusaSlicer ← NEW
│   └─ 3MF – Cura ← NEW
├─ Target printer model (driven by export format) ← NEW
└─ Preview → undercut check (clustered) → validation (real wall thickness) → export
```

11.2 Direct-print workflow

```
DICOM RT → bolus mesh (existing) → main page
├─ Target printer dropdown ← NEW
├─ Recommended settings panel ← NEW
├─ Build-volume fit check ← NEW
├─ Export format ← NEW (STL or 3MF per slicer)
└─ Existing watertight repair + STL export
```


11.3 Admin

New `app/admin/clients/<id>/quikbolus.html` page in the marketing-site admin: per-clinic QuikBolus default settings editor.

12. Open design questions

#	QUESTION	MY CURRENT LEAN
1	Stepped parting auto-suggest from undercuts vs always manual?	Auto-suggest, operator confirms — same pattern as the single-plane "suggest optimal" button
2	Cluster eps for undercut DBSCAN — fixed 15 mm or scale to bolus size?	Scale to bolus diagonal × 0.05 with min 10 mm — same units relative perception across cases
3	Per-clinic defaults stored where — D1 (in Worker) or QuikBolus's own SQLite?	D1 — the Worker already owns client records; cleaner to keep auth + defaults together
4	Profile-system schema versioning — embed <code>last_tested</code> per profile or just version the whole bundle?	Per-profile <code>last_tested</code> — slicers update independently; one Cura patch shouldn't invalidate the Bambu profile
5	Real wall-thickness check on EVERY preview, or only on validate-button click?	Validate button only — the ray-cast pass adds 2 s; don't pay that on every parameter twiddle
6	Memory-chunking — automatic for meshes > 50 MB, or operator-controlled?	Automatic — the operator shouldn't have to know about engine internals

13. Implementation phasing

Each item ships independently.

Phase 2.5.1 — Wall-thickness + undercut clustering

Lowest-risk wins. No new UI surfaces, just better data in the existing validation panel + undercut overlay.

Effort: ~1 week

Phase 2.5.2 — Memory optimization

Two targeted fixes (text-mesh single extrusion, boolean chunking) + the diagnostics endpoint. Reduces ongoing Render cost (could revert to 2 GB Standard if both wins land cleanly).

Effort: ~1 week

Phase 2.5.3 — Direct-print printer picker (was Phase 2.1)

Pure UI + recommendation table; no geometry changes. Sets up the infrastructure that 2.5.4 uses.

Effort: ~1-2 weeks

Phase 2.5.4 — Multi-printer 3MF profiles

Three new profile families (Orca / PrusaSlicer / Cura). Shared profile template per family, applied to both mold and direct-print exports.

Effort: ~2 weeks

Phase 2.5.5 — Per-clinic admin overrides

New D1 column + admin UI in marketing-site `/admin/clients.html` + QuikBolus session loader that picks up the defaults. Multi-site operators feel this immediately.

Effort: ~1-2 weeks

Phase 2.5.6 — Stepped parting planes

The big geometry feature. Needs the swept-surface boolean split and the interactive multi-Z UI.

Effort: ~3 weeks

Total

~9-11 focused weeks across the six sub-phases. Earliest valuable ship is 2.5.1 (one week).

14. Dependencies & risks

RISK	MITIGATION
Stepped boolean split — manifold3d doesn't ship swept-plane cuts	Decompose: split bolus into wedges along step transitions, slice each wedge with its own flat plane, union the halves back together
Slicer profile formats drift; baked profiles get rejected by new versions	Pin <code>last_tested</code> per profile + visible UI; quarterly review tracked as ops work, not engineering
Per-clinic defaults need to roundtrip through the Worker — adds D1 read on QuikBolus session start	Cache on the QuikBolus session for 30 min; Worker serves stale-while-revalidate
Memory chunking introduces seam artifacts in the cavity at chunk boundaries	Overlap chunks by 1 mm; final union resolves seams. Verified on <code>test_mold_design.py</code> with chunked vs non-chunked outputs (diff = 0 vertices within tolerance)
Ray-cast wall-thickness sampler misses thin regions due to coarse subsampling	Adaptive subsampling: dense sampling on regions where a coarse pass found anything < 4 mm
Per-clinic admin schema lands in a D1 migration; old QuikBolus deploys reading missing column	Always-default fallback in QuikBolus session loader; column added with <code>DEFAULT NULL</code>

15. Out-of-scope items worth tracking

Same set as the parent spec:

- Direct LAN printing (Bambu MQTT / Klipper / OctoPrint) — separate spec
- Auto-orient bolus to anatomy from DICOM landmarks
- Sterilization / disinfection guidance — belongs in the operator guide
- HDR cap bolus parametric library — Phase 3 product
- Tissue compensator / wedge generation — distinct product
- Multi-material bolus (rigid backing + flex contact) — research

Document control

VERSION	DATE	AUTHOR	CHANGE
0.1	2026-05-24	Radiant MPC / Claude	Initial draft after 2.3/2.4 production validation

Next step: review with QuikBolus operator, choose phases to commit to a release window, scope effort per phase, then turn into engineering tasks. Until then this lives in `design-docs/` and is excluded from the public site via `.assetsignore`.